

## Product details

### Application of facial mask , breathing filter mask

- 1) Widely used in the daily life,dust free work shop,food factory, to prevent daily life dust ;
- 2) For hospital,clinic,doctor and patients to prevent virus,body fluids and blood splash content.

### Specifications of doctor and nurse face mask

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing	50 pcs/box , 2000 pcs/ctn, or as your requirement
Packing design	all inner box and outer carton printing can be as your design

### Storage of mask

- 1) The product will maintain its properties when it is stored in a dry place at temperatures between 10 ° and 30 ° C.
- 2) Not exposed to direct sunlight.

## Product display



### The regular shipment way



## The happy time with customers



## FAQ

Q: 1.Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003,

we have our own exported rights company in 2013.

Q: 2.What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q:3. What is the number of employees at your company?

A: Our factory is about 380-400 skilled workers, 10 professional QC for quality checking during each process

and 10 professional int'l sales with 3-8 years experience in int'l trade department.

Q:4. What is the percentage of your products distribution for local, and export?

A: At present, 90% for export,10% for local.