

Product details

Usage of patient bed transfer sheet

It's widely used in hospital or clinic, for moving or re-positioning a patient from a bed to gurney.

Non woven sheet with handles cut out at sides.

The feature of patient bed transfer sheet

Fabric	Polypropylene,SMS,PP+PE
Style	Ultrasonic heat sealing,with hole in the length side
Function	Protect staff and patients from cross-infection
Size	195cm*85cm (can be made according to client's demand)
Color	White,blue,green,gray,black....etc
Packaging	Sterile: 1 pcs/paper plastic bag,50 pcs/Ctn Non-sterilize: 1 pcs/bag, 50 pcs/ctn Or as your requirements

Product display



orient Honest
MEDICAL INTERNATIONAL

