

China best manufacturer for earloop face mask , disposable face masks for allergies , medical mask anti virus

Specification of earloop face mask

- 1) Size: 17.5 x 9.5 cm, one size fit the most
- 2) Adaptable nose bar, no stimulus to human skin
- 3) Colors: White/pink/green/blue/yellow, we have 11 colors totally
- 4) High filter efficiency, filter pollen, dust and bacteria, BFE% > 99%
- 5) Resistance to penetration by synthetic blood-120mmHg
- 6) Conform to EN14683 and FDA quality standard

The fabric material of face masks for allergies

- 1) It's non-woven, non-toxic, No irritation, Soft and comfortable ;
- 2) Inner is the PP non-woven with antibacterial effect ;
- 3) Middle is the Efficient filter layer (high-quality Melt blown fabric) ;
- 4) Outer is the Health Grade Print Spun bond Polypropylene .

Packaging of mask anti virus

| Specification | Normal packing | Carton size | 1*40 HQ container can hold (Pcs) |
|---|--------------------------|-------------|----------------------------------|
| 17.5*9.5cm, 3ply, 18+20+20G | 50 pcs/box, 40 boxes/ctn | 52*19*42cm | 3,280,000 Pcs |
| 17.5*9.5cm, 3ply, 25+25+25G | 50 pcs/box, 40 boxes/ctn | 51*37*30cm | 2,400,000 Pcs |
| Remark: all can be special customized according to your actual requirements | | | |

Picture of earloop face mask



orient Honest
MEDICAL INTERNATIONAL

