

Product details

Features of China manufacturer for SMS surgical gowns

1. All-in-one, back opening.
2. The surgical gown is durable, soft, light, comfortable, abrasion and static electricit
3. With perfectly smooth surface, surgical gown is able to expel inorganic liquid and prevent splashing of chemicals. Solid dust cannot attach to the gown as well.
- 4.The surgical gown is specially coated to effectively prevent harmful, superfine dust, high density inorganic acid, alkali and salt solutions.
5. It is proved that the surgical gown provides protection against the attack of blood, virus.

Specification of China manufacturer for SMS surgical gowns

Description	China manufacturer for SMS surgical gowns,the best plant hospital sterile gowns,disposable blue surgery gown supplier
Material	SMS (Polypropylene Spunbond Nonwoven Fabric + Meltblown Nonwoven fabric + Polypropylene Spunbond Nonwoven Fabric)
Size	S(110*130cm), M(115*137cm), L(120*140cm) XL(125*150cm) 2XL(130*150cm) or any other customized sizes
Weight	16-65gsm is available(as your request)
Details	Sterile and non-sterile, standard and reinforced style, standard and ultrasonic, with and without towels. Hook And Loop on collar, ties on waist, knitted cuffs
Color	blue
Packaging	non sterile,10 pcs/bag, 50pcs/ctn or 100pcs/ctn. sterile,1 pcs/sterile bag, 50pcs/ctn or 100 pcs/ctn
Application	Medical & health / Household / Laboratory
Characteristic	Durable, Eco-friendly, Non-toxic, Breathable and Flexible
	Anti-alcohol,Anti-Blood, Anti-Oil, Anti-bacteria,
	Waterproof,Acid proof, Alkali proof, UV-stabilizer, Non-milder
Delivery	One 40HQ container can be finished within 15-25days.
Storage	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Terms of payment	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Note	1.This product can be sterilized by ethylene oxide.

	2.Validity: 3 years
OEM	1.Material or other specifications can be according to customers'requirements. 2.Customized Logo/brand printed. 3.Customized packaging available.
Certificate	CE, ISO13485, FDA

Product display



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

