

Product details

Key features include

- 1) patient transfer sheet is easy-slide design,Ultrasonic heat sealing
- 2) The performance of the fabric is immensely strong,Polypropylene 100gsm up, or PP+PE 100gsm up.
- 3) The transfer slide sheet can hold weight 300 LBS up, facilitates safe and easy transfer of patients.
- 4) Lightweight, latex free,compact,and easy to store
- 5) Single use-protecting staff and patients from cross-infection, cost effective
- 8) The regular size of single-use non woven transfer sheet : 195cm*85cm, or as your requirements.

The specification sheet

Material	PP,PP+PE
Color	White, blue, green....etc
Style	Ultrasonic heat sealing,with hole in the length side
Size	195cm*85cm (can be made according to client's demand)
Application	Hospital,emergency center
MOQ	5,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	10,000pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	According to the quantity,normally 30days after deposit
Packing details	Sterile: 1 pcs/paper plastic bag,50 pcs/Ctn Non-sterilize: 1 pcs/bag, 50 pcs/ctn Or as your requirements


Attention of single-use non woven transfer sheet :

- 1) Store it in dry and cold place ;
- 2) Please stop using if it's dirty ;
- 3) Using in hospital or hotel where people want to save cleaning time.

Product display



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

