

## Product details

### China manufacturer waterproof hospital CPE gown:

China manufacturer waterproof hospital CPE gown is made from hydrophobic CPE material□ Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent. and can be used in hospital, laboratory, food industry, or any where need protection.

### Specification of China manufacturer waterproof hospital CPE gown:

material	CPE
Weight	CPE 30-65 gsm; (comfortable,without glass fibres,latex free.)
Style	Neck: with a hool Waist: ties on Cuff: with thumb loop Opening on back
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Color	blue
packing	non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
packing design	All the inner box and carton printing can be as your design
carton size	Standard packing carton 58X28X48 cm
Applications	hospital clinical medical personnel and patients. used in hospital, laboratory, food industry, and any where need protection

### Product Advantages:

1. This product is Latex-free; abrasion-resistant;completely impervious fluid and alcohol repellent.
2. Disposable isolation gown is designed with a major concern for comfort and performance.
3. Ties on neck and waist,ensuring that gowns remain securely in place to protect clothes.
5. It can reduce the risk of contamination and transmission of infections.


## Product display



orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

**Q: 1. Are you a manufacturer or a trading company, and when was your company established?**

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

**Q: 2. What is the annual turnover for your company?**

A: Normally it is around USD18 million-20 million, it is still in the growing.

**Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?**

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

**Q: 4. Can I get the free samples for quality checking? How long time can I get it?**

A: Yes, our samples can be supplied free within 3-7 days.