

## Product details

### Introduction of green scrubs surgical caps

Green scrubs surgical caps are made of PP nonwoven material, hand made, with ties on the back of the head, strong and adjust ties make the cap very comfortable to wear it. Mainly use for hospital doctor operate room.

### Specifications of green scrubs surgical caps

Material	SBPP/SPP/nonwoven, SMS
Weight	PP10-30g, or SMS35-65g
Style	Elastic (latex free), single and double elastic style
Size	18",19",21",24"28"(as your require)
Color	green, or as your request
Packing	100pcs/bag, 20bags/ctn, 2000pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58x25x27cm
Applications	Hospital clinical medical personnel and patients. Dust-free workshop, Laboratory, food industry, electronic manufactures and so on
MOQ	100,000 PCS

### Indications of green scrubs surgical caps

Sterile green scrubs surgical caps are used for hospital clinical medical personnel and patients;  
Non-sterile green scrubs surgical caps can be used for dust-free workshop, hospital, laboratory, food industry, electronic manufactures and so on.

## Product display



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL SUPPLY



## Our business process



## The regular shipment way



## FAQ

**Q: 1. Are you a manufacturer or a trading company, and when was your company established?**

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

**Q: 2. What is the annual turnover for your company?**

A: Normally it is around USD18 million-20 million, it is still in the growing.

**Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?**

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

**Q: 4. Can I get the free samples for quality checking? How long time can I get it?**

A: Yes, our samples can be supplied free within 3-7 days.

