

## China surgical white ear loop disposable face mask for food service

### Disposable face face mask indications

- 1.Products are used for hospital clinical medical personnel,patients,dust-free workshop, hospital,laboratory, food industry,electronic manufactures and so on.
2. To prevent the spread of body fluids and blood splash content and isolate dust, particle, alcohol, blood, bacteria and virus invading.

### Disposable face mask normal packing

50pcs/box,2000pcs/ctn, carton size:52x38x30cm,  
1\*20ft load:500 cartons,1\*40HQ load:1150 cartons.

### Disposable face mask Specifications

Item	China surgical white ear loop disposable face mask for food service
Material	3ply, pp non woven 20g+filter 20g+pp non wovn 25g
Weight	20+20+25g (comfortable,without glass fibres,latex free.)
Style	Elastic ear-loop (latex free)
Size	17.5x9.5cm(adult)
Color	white
Packing	Non-sterile:50 pcs/box,2000pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	52x38x30cm
Applications	hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on
MOQ	5000pcs
Sample	Can be supplied free for your quality checking within one week
Payment term	T/T,L/C,D/A,D/P/ESCROW
Delivery time	as you quantity(we can finish 1 million pcs per day)
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou