

## Product details

### The feature of personalized face masks

- 1) Soft, comfortable,folds up for easy carrying
- 2) Fluid resistant
- 3) 99% bacterial filtration efficiency
- 4) Standard Ear-Loop

### Our pink surgical mask can achieve following quality standard

- 1) Conform to EN14683: 2005
- 2) Bacterial filtration efficiency (BFE%) > 99%.
- 3) TYPE IIR standard, resistance to penetration by synthetic blood-120mmHg.

### The packaging of safety face mask

Description	Qty/box	Qty/Ctn	Box size: 18*10*7cm
25+25+25g, 3ply	50 Pcs	2000 Pcs	Carton size: 52*38*30cm

## Product display



## The regular shipment way



## Our business process

