

## **Chinese workshop for disposable medical mask , hospital disposal mask , 3ply surgery face mask**

### **Indications of disposable medical mask**

- 1) Sterile face mask is for hospital, clinical, medical personnel and patients;
- 2) Non-sterile face mask can be used for dust-free workshop, food factory and so on.

### **We produce following kinds of hospital disposal mask**

- 1) 3ply or 2 ply non woven face mask
- 2) 2 ply or 1ply paper face mask
- 3) 4 ply active carbon face mask

### **Performance of 3ply surgery face mask**

- 1) To prevent the spread of body fluids and blood splash content.
- 2) Conform to EN14683 and FDA
- 3) Bacterial filtration efficiency (BFE%) > 99%.

### **Packaging of 3ply surgery face mask**

Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
17.5*9.5cm, 3ply, 18+20+20G	50 pcs/box, 40 boxes/ctn	52*19*42cm	3,280,000 Pcs
17.5*9.5cm, 3ply, 25+25+25G	50 pcs/box, 40 boxes/ctn	51*37*30cm	2,400,000 Pcs
Remark: all can be special customized according to your actual requirements			

### **Specifications of disposable medical mask**

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing design	all inner box and outer carton printing can be as your design
MOQ	100,000 pcs

Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day,1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

### Picture of hospital disposal mask





orient Honest  
MEDICAL INTERNATIONAL

