

## Product details

### Specifications of tie on surgical face masks , doctor and nurse mouth covers

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing	50 Pcs/box, 2000 Pcs/Ctn, Or as your requirements
Packing design	all inner box and outer carton printing can be as your design
Applications	hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on

### Storage of breathable face mask with nose clip

- 1) The product will maintain its properties when it is stored in a dry place at temperatures between 10 ° and 30 ° C.
- 2) Not exposed to direct sunlight.

## Product display



orient Honest  
MEDICAL INTERNATIONAL



## The regular shipment way



## Our business process

