

## Product details

### The specifications of tie on mask

Material	PP non woven + high filtration paper + PP non woven
Color	White/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy
Ear loop	Tie on
Nose bar	Plastic or metal
Quality standard	TYPE IIR,EN14683 and FDA

### The feature of sanitary mask

- 1) Three-dimensional design: specified Anti-SARS masks,sterilization,anti-flu,dust
- 2) With nose clip: flexible,fitting fixed
- 3) Fine edge: the edge of full heat sealing technology package
- 4) High elastic cord: hypo-allergenic material

### Our trade terms of hospital masks for sale


- 1) MOQ: 100,000 pcs
- 2) Payment terms: T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
- 3) Production capability: 1\*40 HQ container can be finished in 15 days
- 4) Delivery time: normally, 30 days after deposit
- 5) Sample: sample is free of charge, we can offer samples in 3 days

## Product display



orient Honest  
MEDICAL INTERNATIONAL





## The regular shipment way



## Our business process

