


The Specification of COOL FEELING disposable face mask

Fabric Material	Layer 1: super soft spunlace Layer 2: high filtration paper Layer3: ES material
Style	Ear loop
Dimensions size	17.5*9.5 cm
Fabric weight	25gsm + 25gsm + 35gsm, could be as your request
Color	white, baby blue, baby pink, gray
Certificate	CE, ISO13485, FDA
MOQ:	100000pcs
Packing design	can be as your design
Packing	50pcs/box,2000pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	7-30 days
Applications	hospital, laboratory, clinic, food industry, dust-free workshop


The picture of COOL FEELING disposable face mask



Precision welding
resistant to pull-off



Antiparticulate
microorganism




Anti-virus, anti-bacteria



orient Honest
MEDICAL INTERNATIONAL

0.32 cool sense of value



test item			
instant cool value	0.2-0.3	0.3-0.4	>0.4
somatosensory	not obvious	cool without ice	Have a cold feeling

* Data from the laboratory, different colors, fabric batches will have a certain impact on the value