



|                       |   |
|-----------------------|---|
| Product Name          | 3D Fish Shaped KF94 Face Mask                             |
| BFE,PFE               | BFE≥95%   |
| Material              | PP+ Melt Blown+Melt Blown+ PP                             |
| Color                 | Multiple colors, or printed, or as per request            |
| Style                 | Ear-loop  |
| Weight                | 50g+25g+25g+30g, or as your request                       |
| Size                  | 20x8cm for Adults, contact for kids' size                 |
| Packing               | 50pcs/bag, 40bags/carton, 2000pcs/carton                  |
| Delivery Time         | 7-15 days   |
| Terms of Payment      | T/T, L/C, X-transfer, PayPal, etc                         |
| Service               | OEM, Samples provided, Provision of Product Specification |
| Testing Certification | GB2626-2019   |

Samples can be provided as requested. Welcome your enquiry.

- \* 4 layers providing better protection
- \* 3D style creating more room for cooler and fresher breathing
- \* Appealing looking, lightweight, closer touch of mask edge to the skin
- \* Sensitive skin-friendly
- \* Less stress ear-loop on the ear











Available in different weights, sizes, colors and packaging per requested