

## Application of disposable ear-loop face mask

- 1) For hospital clinical medical and patients,dust free work shop,food factory.
- 2) To prevent of spread of body fluids and blood splash content.

## Features of disposable hospital masks

- 1) Disposable ear-loop face mask is made of 3-layers,inside and outside layers are PP non woven,middle layer is high standard N95/N98 filter paper.
- 2) The nose clip can be adjusted comfortable according to the different face shapes;  
The ear loop is firm and not easy to fall;
- 3) Non-irritating and soft,Avoiding bacteria,dust,foreign body into pneogaster,
- 4) Remove harmful gas, peculiar smell, stench.
- 5) BFE is over 99%, CE/ISO13485, type II /Type IIR standard.

## Packaging of disposable mask

Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
17.5*9.5cm,3ply, 18+20+20G	50 pcs/box,40 boxes/ctn	51*19*41cm	3,450,000 Pcs
17.5*9.5cm,3ply, 25+25+25G	50 pcs/box,40 boxes/ctn	51*37*30cm	2,400,000 Pcs
Remark: all can be special customized according to your actual requirements			

## Specifications of disposable mask

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ea rloop
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue /burgundy/black
Packing	50 pcs/box,2000 pcs/ctn, or as your requirement
Packing design	all inner box and outer carton printing can be as your design
Applications	hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on

## Why do customers choose us ?

- 1) We are the direct disposable ear-loop face mask supplier China, factory price.

- 2) Stable quality, CE,ISO certification.
- 3) Production capacity: 1000,000pcs/day, 1\*40HQ can be finished within 10-15 days;
- 4) OEM packing design, service in 24 hours.

#### Picture of disposable hospital masks





orient Honest  
MEDICAL INTERNATIONAL

