

Product details

Parameters of medical child face mask

- 1) Material: PP non woven + high filtration paper + PP non woven, very soft, no irritation.
- 2) Color: White
- 3) Ear loop: elastic flat ear loop, 6mm width
- 4) Size: 125×80mm
- 5) Nose bar: no nose bar
- 6) BFE(Bacterial filtration efficiency) > 99%

Packaging of nonwoven medical child face mask

		Carton size	1*40 HQ container can hold (Pcs)
Normal packaging	50 Pcs/box, 2000 Pcs/Ctn	52*29*34cm	2,620,000 Pcs
Single packaging	1 Pcs/1 color printed sealed plastic bag, 50 Pcs/box, 2000 Pcs/Ctn	52*30*42cm	2,060,000 Pcs

Pictures of disposable medical child face mask



orient Honest
MEDICAL INTERNATIONAL

