

Nonwoven tie on surgical mask , 3 plys surgical face mask , factory for hospital masks

Properties of tie on surgical mask

- 1) It's very low resistance to breathing, free from foreign smell, without any stimulating feeling;
- 2) Nose bar can be adaptable, comfortable, without glass fibers, latex free;
- 3) BFE (Bacteria filtration efficiency) $\geq 99\%$
- 4) Use ultrasonic welding technology, comfortable wearing.

Packaging of tie on surgical mask

50 pcs/box, 2000 pcs/ctn

Box size: 18*10*10cm

Ctn size: 52*38*42cm

Our surgical face mask can achieve following quality standard

- 1) Conform to EN14683: 2005
- 2) Bacterial filtration efficiency (BFE%) $> 99\%$.
- 3) TYPE IIR standard, resistance to penetration by synthetic blood-120mmHg.

Specifications of hospital masks

Material	PP non woven + high filtration paper + PP non woven
Color	White/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Ear loop	Tie on
Nose bar	Plastic or metal
Quality standard	TYPE IIR, EN14683 and FDA
MOQ	100,000 pcs
Price:	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	1,000,000 pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Packing details	50 pcs/box, 2000 pcs/ctn, or as your requirement

Picture of surgical face mask



orient Honest
MEDICAL INTERNATIONAL

