

## Product details

### Indications of patient transfer sheet

- 1) The patient transfer sheet is used to re-position the patients.
- 2) Recommend single patient use.
- 3) Heavy weight non woven fabric with treated inner surface that assist caregivers with re-position tasks.
- 4) Glides easily in all directions and significantly reduces friction and force required to move patients.

### The feature of patient transfer sheet

Fabric	Polypropylene,SMS,PP+PE
Style	Ultrasonic heat sealing,with hole in the length side
Function	Protect staff and patients from cross-infection
Size	195cm*85cm (can be made according to client's demand)
Color	White,blue,green,gray,black....etc

### Packing details of patient transfer sheet

Sterile: 1 pcs/paper plastic bag,50 pcs/Ctn

Non-sterilize: 1 pcs/bag, 50 pcs/ctn

Or as your requirements

## Product display



orient Honest  
MEDICAL INTERNATIONAL

