

## Product details

### Description of anti flu face masks , medical facemasks

- 1) Top safe,3-Ply, fluid resistant layer with unique single pleated design
- 2) Soft and elastic ear-loop, round or flat
- 3) Adjustable nose piece, metal or plastic material
- 4) Fiberglass-free,latex-free, BFE > 99.8%, PFE > 99.3%
- 5) Packaging: 50 Pcs/box, 2000 Pcs/Ctn

### Specifications of earloop procedure masks

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing design	all inner box and outer carton printing can be as your design
Applications	hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on

## Product display



### The regular shipment way



## Our business process

