

Product details

The usage of Polypropylene medical masks

- 1) The mask is widely used in hospital, clinical, medical personnel and patients.
- 2) It offer first-rate protection against the dangerous air found in all industrial work situations such as painting,
lab work involving chemical solvents, spraying chemical fertilizer or pesticides, odor guard when doing waste management,
working with oils and gas fumes...etc.

Specification of antiviral facemask

- 1) Size: 17.5 x 9.5 cm, one size fit the most
- 2) Adaptable nose bar, no stimulus to human skin
- 3) Colors: White/pink/green/blue/yellow, we have 11 colors totally
- 4) High filter efficiency, filter pollen, dust and bacteria, BFE% > 99%
- 5) Resistance to penetration by synthetic blood-120mmHg
- 6) Conform to EN14683 and FDA quality standard

The packaging of BFE99 masks

Q
b
7
6
b
b
p
t
1
8
n
1
*
7
c
m
C
2
B
t
2
b
4
2
5
9
;
3
p
2
y
3
8
*
3
c
m

Product display



The regular shipment way



Our business process

