

Product details

Advantage of custom hospital gowns light blue

- 1.PP material: Lint-free, waterproof, good tensile strength, comfortable.
- 2.Light-weight, Non-irritating, antibacterial property ,good tensile strength, Superior breathability and anti-static, impact resistance, soft, comfortable.

Specification of custom hospital gowns light blue

Name	custom hospital gowns light blue,dental disposable isolation gown,disposable dental gowns manufacturer
Material	Polypropylene,PP Nonwoven Fabric
Style	Non-sterile, standard style, Tie on collar, ties on waist, knit cuffs
Size	S(110*130cm), M(115*137cm), L(120*140cm), XL(125*150cm), 2XL(130*150cm) or any other customized sizes
Grammage	25g
Color	light blue
Packaging	non sterile,10 pcs/bag, 50pcs/ctn or 100pcs/ctn,
Container loading	20'FCL:680CTNS,40'HQ:1700CTNS
MOQ	5000 PCS
Delivery	One 40HQ container can be finished within 15-25days
Store	Store in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

Product display




orient Honest
MEDICAL INTERNATIONAL





Our business process



The regular shipment way



FAQ

Q: 1.Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2.What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE,ISO13485,FDA. REACH and NELSON report.

Q: 4.Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

The happy time with customers

