

## Product details

### Introduction of disposable SMS dental surgical lab coat

1. Material: SMS
2. Color: green, etc.
3. Grammage: 50gsm
4. Style:knitted Collar and cuffs, with three pockets
5. Size:S-2XL, or custom,

### Specifications of disposable SMS dental surgical lab coat

Name	disposable SMS dental surgical lab coat,custom wholesale surgical lab coats,hostpital disposable lab coats
Material grammage	SMS50g
Style	with knit collar and cuffs, with three pockets
Size	S M L XL 2XL
Color	green, or blue, yellow, pink, red, grey, black etc.
Packing	Non sterile:5pcs/bag,10 bags/ctnor ; Sterile:1pcs/bag, 50pcs/ctn
Packing design	All inner box and carton printing can be as your design
MOQ	5,000 pcs
Sample	Can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

## Product display



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL





## Our business process



## The regular shipment way



## FAQ

**Q: 1. Are you a manufacturer or a trading company, and when was your company established?**

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

**Q: 2. What is the annual turnover for your company?**

A: Normally it is around USD18 million-20 million, it is still in the growing.

**Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?**

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

**Q: 4. Can I get the free samples for quality checking? How long time can I get it?**

A: Yes, our samples can be supplied free within 3-7 days.

## The happy time with customers

