

## Pattern of disposable hospital isolation gowns

1. Long sleeves
2. Elastic cuffs
3. Ties on neck and waist
4. With sterilized and non-sterilized packing
5. Different specifications are available

## Indications of disposable hospital isolation gowns

1. Prevent the spread of body fluids and blood splash content and isolate dust, particle, alcohol, blood, bacteria and virus invading.
2. Perfect for operations with large amounts of fluids. For preventing a transfer of infectious from or staff to the wound or vice versa by direct contact.

## Advantages of disposable hospital isolation gowns

1. Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent.
2. Providing reliable barriers to water steam or blood and other fluids.
3. Can prevent pollution and the surface of this clothes without stains.
4. Isolation gowns are fluid resistant and offer economical, comfortable and reliable protection.
5. These isolation gowns have elastic cuffs with sewn seams which provide superior strength.
6. It can reduce the risk of contamination and transmission of infections.

## Specifications of disposable hospital isolation gowns

Material	PP non woven
Style	Tie/velcro on collar, tie on waist, elastic/knitted cuff
Size	115*137cm,120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm,etc
Color	White, green, blue, yellow, pink, red, grey, black etc.
Carton size	58X28X48 cm
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Weight	PP 20-65gsm; PP+PE 40-65 gsm; SMS 35-65 gsm (comfortable, without glass fibres, latex free.)
Applications	Hospital clinical medical personnel and patients. dust-free workshop, laboratory, food industry, electronic manufactures and so on

## Pictures of disposable hospital isolation gowns





orient Honest  
MEDICAL INTERNATIONAL

