

Product details

Features of disposable surgical gown reinforced

1. The surgical gown is durable, soft, light, comfortable, abrasion and static electricit
2. With three defense function:anti blood, anti alcohol, anti static.
- 3.The surgical gown is specially coated to effectively prevent harmful, superfine dust, high density inorganic acid, alkali and salt solutions.
- 4.It is proved that the surgical gown provides protection against the attack of blood, virus.

The packaging of disposable surgical gown reinforced:

Sterile: 1 pcs/paper plastic bag,50 pcs/Ctn

Non-sterilize: 5 pcs/bag, 50 pcs/ctn

Or as your requirements.

Specifications of disposable surgical gown reinforced

Name:	disposable surgical gown reinforced,SMS surgical gown with reinforced layer,China surgical gown with ties for sale
Material:	SMS
Grammage:	50 gsm
Closure:	with belt closure, or with Velcro closing
Application:	for medical, surgical, pharmaceutical, laboratorial, industrial sector
Size:	S/M/L/XL/XXL/XXXL,
Color:	Green,Blue, etc, optional
Packing:	Normal packing:1pcs/polybag,50pcs/carton
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Certification:	CE, ISO13485,FDA

Product display



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

