

### Description of medical surgical scrub hats:

Product Name:	doctor surgeon cap manufacturer China,hospital surgical bouffant caps,medical surgical scrub hats	Material:	PP Nonwoven
Color:	white (other Colors Can Be Available)	Style:	single elastic
Size:	19"	Grammage:	PP10g

### Advantages of medical surgical scrub hats:

- 1.Pleated elastic band around, ample room for long hair. It can prevent hairs falling to avoid possible troubles.
- 2.Used in hospitals and clean room.
- 3.Breathable, non-absorbent, durability
- 4.PP nonwoven, economical, elegant and good looking.

### Indications of medical surgical scrub hats:

The medical surgical scrub hats are used for hospital clinical medical personnel and patients. And the non-sterile disposable caps can be used for dust-free workshop,, hospital, laboratory, food industry, electronic manufactures and so on.

### Specifications of doctor surgeon cap manufacturer China:

Material	PP non woven
Grammage	10g
Style	single elastic
Size	19", etc.
Color	white, (other colors can be available)
packing	100pcs/bag, 20bags/ctn, 2000pcs/ctn
packing design	All inner box and carton printing can be as your design

carton size	58x25x27cm
Applications	For hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on
MOQ	100,000 PCS

**Detail pictures of doctor surgeon cap manufacturer China:**



