

## Product details

### Introduction of exporter for disposable CPE protective gown:

Exporter for disposable CPE protective gown is made from hydrophobic CPE material, waterproof, Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent. and can be used in hospital, laboratory, food industry, etc.

### Specification of vendor for disposable CPE isolation gown:

material	CPE
Weight	CPE 30-65 gsm; (comfortable, without glass fibres, latex free.)
Style	Neck: with a hood Waist: ties on Cuff: with thumb loop Opening on back
Size	120x130cm, 130x150cm, 150x160cm, 150x190cm, 170x170cm
Color	milky white
packing	non sterile, 10 pcs/bag, 10 bags/ctn. sterile, 1 pcs/sterile bag, 100 pcs/ctn
packing design	All the inner box and carton printing can be as your design
carton size	Standard packing carton 58X28X48 cm
Applications	hospital clinical medical personnel and patients. used in hospital, laboratory, food industry, and any where need protection

## Product display



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

**Q: 1. Are you a manufacturer or a trading company, and when was your company established?**

**A:** Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

**Q: 2. What is the annual turnover for your company?**

**A:** Normally it is around USD18 million-20 million, it is still in the growing.

**Q: 3. What is the number of employees at your company?**

**A:** Our factory is about 380-400 skilled workers, 10 professional QC for quality checking during each process and 10 professional int'l sales with 3-8 years experience in int'l trade department.

**Q: 4. What is the percentage of your products distribution for local, and export?**

**A:** At present, 90% for export, 10% for local.

## The happy time with customers

