

## Product details

### Features of SMS surgical gown:

- 1.SMS material: Lint-free, waterproof, good tensile strength, comfortable.
- 2.Light-weight, Non-irritating, antibacterial property ,good tensile strength, Superior breathability and anti-static, impact resistance, soft, comfortable.
- 3.Provide an excellent barrier to potentially infectious blood and other bodily fluids. against splatter of chemical material and solid dust.
- 4.Disposable SMS surgical gown is the most conventional medical apparel, can reach the health protection standards in the general case.
- 5.SMS : “Three defenses” (waterproof, anti-alcohol and anti-blood spatter) to comply with domestic and international medical standards.

### Details Information of SMS reinforced surgical gown:


Name:	green SMS reinforced surgical gown plant, hospital sterile operating gown, sterile packing reinforced surgical gown
Material:	SMS
Grammage:	50 gsm
Closure:	with belt closure, or with Velcro closing
Application:	for medical, surgical, pharmaceutical, laboratorial, industrial sector
Size:	S/M/L/XL/XXL/XXXL
Color:	Green,Blue, etc, optional
Packing:	Normal packing:1pcs/polybag,50pcs/carton
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Certification:	CE, ISO13485,FDA

## Size Chart:

### SURGICAL / ISOLATION GOWN SIZE CHART (Customized size chart is welcome)

CM	S	M	L	XL	XXL	XXXL
HEIGHT	115	115	120	125	130	135
CHEST LENGTH	127	137	140	145	150	155

## SMS diagram:



## Product display



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

- 1, Free sample.
- 2, Small order acceptable;
- 3, L/C, T/T, D/P terms available;
4. FOB, CNF, CIF, EXW terms available;
- 5, Competitive and stable price because we are factory.
- 6, Qualified products. 3rd part product inspection is acceptable.
- 7, OEM service. More options on material, pattern, size and packing;
- 8, Timely delivery. The best lead time is 25-30 day for a 40HQ container.
- 9, Professional export service. Documentation, clearance, sea transportation, all in one.
- 10, Good after sale service. We take care complaints as we would like to be your long-term supplier!