

Product details

Features of medical doctor surgery gowns

Medical doctor surgery gowns for operating room , disposable lab coats are made from hydrophobic polypropylene material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency.

Material grams:25-65gsm.

Quality Control

- 1.We will not start to make the products until you confirm the sample.
- 2.The packing will be water-proof, damp-proof and sealed .

Sample Policy

- 1.Sample according to your design drawing. Sampling time:3-5days .
2. Design your own sample, then remake it.
3. Existing Samples Sampling time:1day
4. Make mass production after you confirm the sample.

Specification

| | |
|----------------|--|
| Material | PP+PE |
| Weight | 35gsm ,etc |
| Color | green, etc. |
| Size | 120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm |
| Style | Tie on / Hook And Loop on the neck, elastic / knitted cuff, standard or reinforced |
| Packing | Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn |
| Packing design | All inner box and carton printing can be as your design |
| Carton size | 58X28X48 cm |

Product display




orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. What is the number of employees at your company?

A: Our factory is about 380-400 skilled workers, 10 professional QC for quality checking during each process and 10 professional int'l sales with 3-8 years experience in int'l trade department.

Q: 4. What is the percentage of your products distribution for local, and export?

A: At present, 90% for export, 10% for local.