

## Product details

### Features of medical doctor surgery gowns

Medical doctor surgery gowns for operating room , disposable lab coats are made from hydrophobic polypropylene material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency.

Material grams:25-65gsm.

### Quality Control

- 1.We will not start to make the products until you confirm the sample.
- 2.The packing will be water-proof, damp-proof and sealed .

### Sample Policy

- 1.Sample according to your design drawing. Sampling time:3-5days .
2. Design your own sample, then remake it.
3. Existing Samples Sampling time:1day
4. Make mass production after you confirm the sample.

### Specification

Material	PP+PE
Weight	35gsm ,etc
Color	green, etc.
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Style	Tie on / Hook And Loop on the neck, elastic / knitted cuff, standard or reinforced
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58X28X48 cm

## Product display



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INSTRUMENTS



## Our business process



## The regular shipment way



## FAQ

**Q: 1. Are you a manufacturer or a trading company, and when was your company established?**

**A:** Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

**Q: 2. What is the annual turnover for your company?**

**A:** Normally it is around USD18 million-20 million, it is still in the growing.

**Q: 3. What is the number of employees at your company?**

**A:** Our factory is about 380-400 skilled workers, 10 professional QC for quality checking during each process and 10 professional int'l sales with 3-8 years experience in int'l trade department.

**Q: 4. What is the percentage of your products distribution for local, and export?**

**A:** At present, 90% for export, 10% for local.