

## Introduction of disposable clip cap:

Brand Name: ORIENT HONEST

Model Number: OH-1001, OH1002, OH1003

Certification: CE / ISO13485 / FDA

Place of Origin: China

Material : PP nonwoven

Color: white,green,blue,pink,yellow,pink,orange,purple,red (all colors can be available)

Style: elastic clip caps

Weight: 10g-30g

Size: 18''□19''□21''□etc

Packing details: 100 pcs/bag,20 pcs/ctn

## Specifications of disposable clip cap:

<b>Product:</b>	medical consumables high quality colorful non-woven single or double elastic 19" 20" 21" custom size disposable clip cap
-----------------	--

**Material:**

Nonwoven pp with elastic

**Grammage:**

10g~16g

**Style:**

Single elastic or double elastic

**Size:**

19", or 18",20",21", etc.

**Color:**

White, Green, yellow, Blue, etc.

**MOQ**

100,000pcs

**Packaging:**

100pcs/bag ; 2000pcs/ctn

**Application:**

for household, medical, surgical, pharmaceutical, laboratorial, industrial sectors, painting, cleaning room, etc.

**Storage:**

Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

**Terms of payment:**

T/T, L/C, Paypal,D/P,D/A,etc.

**Other service:**


OEM and Customized Logo/brand

**Certification:** CE, ISO13485,FDA

**Features of disposable clip cap:**

1. Breathable, non-absorbent, durability
2. With single or double elastic available, in pleated or clip shape with various colours available.
3. Bacteria filtration efficiency (BFE) has reached 99.1%. This rate is much higher than normal standard of 97%.
4. Free samples and customization available
5. Clip-shaped cap before unfolded
6. We could guarantee best quality, fast delivery, development new material.

**Pictures of disposable clip cap:**





orient Honest  
MEDICAL INTERNATIONAL

