

## Introduction of disposable clip cap:

Brand Name: ORIENT HONEST

Model Number: OH-1001, OH1002, OH1003

Certification: CE / ISO13485 / FDA

Place of Origin: China

Material : PP nonwoven

Color: white,green,blue,pink,yellow,pink,orange,purple,red (all colors can be available)

Style: elastic clip caps

Weight: 10g-30g

Size: 18''□19''□21''□etc

Packing details: 100 pcs/bag,20 pcs/ctn

## Specifications of disposable clip cap:

### Product:

medical consumables high quality colorful non-woven single or double elastic 19" 20" 21" custom size disposable clip cap

### Material:

Nonwoven pp with elastic

### Grammage:

10g~16g

### Style:

Single elastic or double elastic

### Size:

19", or 18'',20'',21'', etc.

### Color:

White, Green, yellow, Blue, etc.

### MOQ

100,000pcs

### Packaging:

100pcs/bag ; 2000pcs/ctn

### Application:

for household, medical, surgical, pharmaceutical, laboratorial, industrial sectors, painting, cleaning room, etc.

### Storage:

Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

### Terms of payment:

T/T, L/C, Paypal,D/P,D/A,etc.

### Other service:

OEM and Customized Logo/brand

**Certification:** CE, ISO13485,FDA

### **Features of disposable clip cap:**

1. Breathable, non-absorbent, durability
2. With single or double elastic available, in pleated or clip shape with various colours available.
3. Bacteria filtration efficiency (BFE) has reached 99.1%. This rate is much higher than normal standard of 97%.
4. Free samples and customization available
5. Clip-shaped cap before unfolded
6. We could guarantee best quality, fast delivery, development new material.

### **Pictures of disposable clip cap:**



