

**Features of disposable surgery gown for wholesales:**

1. Breathable, softable and strong adsorption ability.
2. Style: with ties or Hook And Loop at neck, and ties on the waist, elastic or knit cuff.
3. Standard size , grammage, color and style are available.
4. Characteristics: environmentally-friendly, good tensile strength, feeling soft, non-toxic, air-permeable. Three defenses: anti static, anti blood, anti alcohol.

**Application of disposable surgery gown for wholesales:**

1. Can be used in hospital, laboratory, food industry, anywhere which need protection, Prevent and isolate dust, particle, alcohol, blood, bacteria and virus invading.
2. Can be used as a visit coat . Prevent the cross infection in medical treatment at sanitation field.


**Specification of disposable surgery gown for wholesales:**

Material	SMS50g
Style	Tie/Hook And Loop on collar, tie on waist, elastic/knitted cuff
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Color	blue, etc.
Carton size	58X28X48 cm
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Grams	SMS 35-65 gsm (comfortable, without glass fibres, latex free.)
Applications	Hospital clinical medical personnel and patients. dust-free workshop, laboratory, food industry, electronic manufactures and so on

**Picture of disposable surgery gown for wholesales:**



orient Honest  
Klasifikasi: 100% Polyester







### **Our advantages:**

- 1. 260 workers, 25000 ㎡ dust free workshop,
- 2. 12 years experience, ISO, CE & FDA certificated,
- 3. One hour reach Wuhan port, nearby Shanghai port,
- 4. Free samples and quick response,

Welcome to email me or call me directly, let's talk about more.

### **Contact information:**

**VIP service**

**Remey Chi**

**Mobile: 0086-15623307739**

**Tel:0086-27 8782 9195**

**E-mail:sales04@orienthonest.com**

**Skype:orienthonest04**