

Features of disposable surgery gown for wholesales:

1. Breathable, softable and strong adsorption ability.
2. Style: with ties or Hook And Loop at neck, and ties on the waist, elastic or knit cuff.
3. Standard size , grammage, color and style are available.
4. Characteristics: environmentally-friendly, good tensile strength, feeling soft, non-toxic, air-permeable. Three defenses: anti static, anti blood, anti alcohol.

Application of disposable surgery gown for wholesales:

1. Can be used in hospital, laboratory, food industry, anywhere which need protection, Prevent and isolate dust, particle, alcohol, blood, bacteria and virus invading.
2. Can be used as a visit coat .Prevent te cross infection in medical treatment at sanitation field.

Specification of disposable surgery gown for wholesales:

Material	SMS50g
Style	Tie/Hook And Loop on collar, tie on waist, elastic/knitted cuff
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Color	blue, etc.
Carton size	58X28X48 cm
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Grams	SMS 35-65 gsm (comfortable, without glass fibres, latex free.)
Applications	Hospital clinical medical personnel and patients. dust-free workshop, laboratory, food industry, electronic manufactures and so on

Picture of disposable surgery gown for wholesales:



orient Honest
Klasifikasi: 100% Polyester





Our advantages:

1. 260 workers, 25000 ㎡ dust free workshop,
2. 12 years experience, ISO, CE & FDA certificated,
3. One hour reach Wuhan port, nearby Shanghai port,
4. Free samples and quick response,

Welcome to email me or call me directly, let's talk about more.

Contact information:

VIP service

Remey Chi

Mobile: 0086-15623307739

Tel:0086-27 8782 9195

E-mail:sales04@orienthonest.com

Skype:orienthonest04