

## Introduction of surgical caps:

Surgical caps made from SBPP nonwoven material, In the standard of CE and ISO13485.

Medical mob caps also called disposable clip or strip cap, are made of PP nonwoven material and pleated elastic band round, this design has ample room for long hair. It can prevent hairs falling then to avoid possible troubles.

## Advantages of surgical caps:

1. Different sizes and colors available.
2. Breathable, non-absorbent, durability
3. With single or double elastic available, in pleated or clip shape with various colours available.
4. Bacteria filtration efficiency (BFE) has reached 99.1%. This rate is much higher than normal standard of 97%.
5. Free samples and customization available
6. Clip-shaped cap before unfolded.

## Specifications of surgical caps :

Material	PP non woven
Weight	10-30g
Style	single elastic, or double elastic style
Size	18", or 19",21",24"28"
Color	green,or white,blue,pink,yellow,pink,orange,purple,red(all colors can be available)
packing	100pcs/bag, 20bags/ctn, 2000pcs/ctn
packing design	All inner box and carton printing can be as your design
carton size	58x25x27cm
Applications	hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on

MOQ

100,000 PCS

**Details picture of surgical caps:**



