

Introduction of reinforced sterile hospital surgery gown:

1. Prevent the spread of body fluids and blood splash content and isolate dust, particle, alcohol, blood, bacteria and virus invading.
2. Perfect for operations with large amounts of fluids. For preventing a transfer of infectious from or staff to the wound or vice versa by direct contact.
3. Long style, with good protective functions.


Specification of reinforced sterile hospital surgery gown

Name	the best Chinese factory for disposable reinforced sterile hospital surgery gown
Material	SMS45g
Style	sterile, standard style, Tie on collar, ties on waist, elastic cuffs
Size	S(110*130cm), M(115*137cm), L(120*140cm), XL(125*150cm), 2XL(130*150cm) or any other customized sizes
Grammage	45g
Color	green.
Packaging	sterile, 1pcs/bag, 50pcs/ctn
Container loading	20'FCL:680CTNS, 40'HQ:1700CTNS
MOQ	5000 PCS
Delivery	One 40HQ container can be finished within 15-25days
Store	Store in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

Pictures of reinforced sterile hospital surgery gown:



orient Honest
MEDICAL INTERNATIONAL





Our advantages:

- 1. 260 workers, 25000 ㎡ dust free workshop,
- 2. 12 years experience, ISO, CE & FDA certificated,
- 3. One hour reach Wuhan port,
- 4. Free samples and quick response,

Welcome email me or call me directly, let's talk about more.

Contact information:

VIP service

Remey Chi

Mobile: 0086-15623307739

Tel:0086-27 8782 9195

E-mail:sales04@orienthonest.com

Skype:orienthonest04