

Product details

Features of disposable sms surgeon gowns :

It is made from hydrophobic SMS material ,Latex-free; abrasion-resistant; low lint;
With a high level of fluid repellency; a good barrier to blood, body fluids and pathogens.
With standard and reinforced style .Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.
Available color: white, blue.

Specification of disposable sms surgeon gowns :

| | |
|----------------|--|
| Name | the best disposable hospital gowns supplier,disposable sms surgeon gowns ,disposable surgical clothes exporter |
| Material gram | SMS 40 gsm (comfortable, without glass fibers, latex free.) |
| Style | Tie or Hook And Loop on collar, tie on waist, knitted cuffs |
| Size | 115*137cm, 120x130cm,130x150cm,150x160cm , 170x170cm or as your requirement |
| Color | white, green, blue |
| packing | Sterile: 1pcs/sterile bag, 50pcs/ctn; Non sterile: 5pcs/bag, 50pcs/ctn |
| packing design | All inner bag and carton printing can be as your design |
| Certifications | CE / ISO13485 / FDA |
| Sample | can be free supplied for your checking within 2 days |
| Applications | Mainly used in hospital, clinical medical personnel and patients. |

Indications:

1. To prevent the spread of body fluids and blood splash content and isolate dust, particle, alcohol, blood, bacteria and virus invading.
2. Perfect for operations with large amounts of fluids. For preventing a transfer of infectious from OR staff to the wound or vice versa by direct contact.

Our Advantage:

1. Samples and customization available.
2. Professional QC to control the product quality. All the products comply with international quality standard such as CE,ISO13485.
3. Have enough experience and gained a high reputation from our customers all the world. Such as America, South America, Europe, Southeast Asia, The Middle East etc.
4. Quick payment, prompt delivery and competitive price.
5. Stable quality and best after-sale service.
6. OEM service is available.

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

The happy time with customers

