

Product details

PP / PE waterproof medical surgical gown

PP/ PE disposable isolation gowns is made from PP+PE material and designed for optimal protection during surgical procedures of moderate duration involving moderate amounts of body fluids and moderate risk of infection.

Specification of waterproof medical surgical gown

Material	PP+PE
Weight	35gsm ,etc
Color	light blue, etc.
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Style	Tie on / Hook And Loop on the neck, elastic / knitted cuff, standard or reinforced
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58X28X48 cm

Product Advantages:

1. Use Eco-friendly material Safety for wearing,no stimulus to skin,light and breathable for body..
2. Soft hand feeling and comfortable. Latex-free; abrasion-resistant; could completely impervious fluid and alcohol repellent.
3. Also can prevent pollution and the surface of this clothes without stains.

Indications of waterproof medical surgical gown

1. Perfect for operations with large amounts of fluids. For preventing a transfer of infectious from OR staff to the wound or vice versa by direct contact.
2. To prevent the spread of body fluids and blood splash content and isolate dust, particle, alcohol, blood, bacteria and virus invading.

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

The happy time with customers

