

Product details

Material features

PP+PE fabric hospital gown for operating room , disposable lab coats are made from hydrophobic polypropylene material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency.

Material weight:25-65gsm.

Quality Control

- 1.We will not start to make the products until you confirm the sample.
- 2.The packing will be water-proof, damp-proof and sealed .

Sample Policy

- 1.Sample according to your design drawing. Sampling time:3-5days .
2. Design your own sample, then remake it.
3. Existing Samples Sampling time:1day
4. Make mass production after you confirm the sample.

Specification of vendor for PP+PE fabric hospital gown

Material	PP+PE
Weight	25gsm ,etc
Color	blue, etc.
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Style	Tie on / Hook And Loop on the neck, elastic / knitted cuff, standard or reinforced
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58X28X48 cm

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1.Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2.What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q:3. What is the number of employees at your company?

A: Our factory is about 380-400 skilled workers, 10 professional QC for quality checking during each process and 10 professional int'l sales with 3-8 years experience in int'l trade department.

Q:4. What is the percentage of your products distribution for local, and export?

A: At present, 90% for export,10% for local.

The happy time with customers

