

## The design details of waterproof disposable white lab coats:

**Collar:** knitted collar or turn down collar (single/double collar )

**Front:** with plastic bottom, Velcro or zipper in front, With / without pocket

**Cuffs:** open, elastic cuffs or knitted cuffs

**Back situation:** The end of the back can be open with about 20cm in order to make the user more comfortable

## Application of waterproof disposable white lab coats:

1. Suitable for isolation and protective use in household, hospital, garden, Painting, food industry and other situations
2. Light weight, cool and breathable; Offer light fluid protection from sprays or mists, as well as dust and dirt.
3. Can be used in hospital, laboratory, food industry, any where need protection.  
Can be used as a visit coat.
4. Sanitation and quality in accordance with the CE, ISO13485 standards.

## Size Chart:

Size	S	M	L	XL	2XL	3XL
Height (CM)	103	111	115	120	124	128
chest circumference(CM)	116	126	132	34	144	154
Sleeve length(CM)	57	57	57	58	58	58

## Specifications of waterproof disposable white lab coats:

Material	Microporous film
Grammage	50g, Could as your request
Style	With/without collar, Snaps or Velcros in front, elastic or knitted cuff, with or without cuff
Size	S-6XL
Color	white, green, blue, yellow, pink, red, gray, black etc.
Packing	5pcs/box, 50 bags/ctn, or as your requirement

Packing design	all inner bag and carton printing can be as your design
Applications	dustfree workshop, hospital, laboratory, food industry, electronic manufactures and so on

### Picture of waterproof disposable white lab coats





orient Honest  
MEDICAL INTERNATIONAL

