

Product details

Product description:

Material: PP (2 layers) & Top V Interlocked

1. Elastic waistband (2 layers) & Elastic
2. Elastic cuffs, 5cm wide
3. No seaming
4. Material: PP
5. Non-woven
6. Disposable
7. Dustless

Product specifications:

Material	workwear, workwear
Size	S-5XL
Weight	PP50g
MOQ	50 suits
MOQ	5,000 suits
Color	White
Features	1. With elastic waistband (2 layers) & Elastic 2. Elastic cuffs, 5cm wide 3. No seaming
Packaging	1suit / polybag, 50suits /
Lead time	1 * 40HQ 15-20days
Payment	T / T, L / C, D / A, D / P, Western Union, Paypal
Certification	ISO13485 CE, FDA, Type5 & 6
Usage	Biohazard

Product display



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are your products made in China?

A: Yes, our products are made in China. We have been producing since 2003 and have a good reputation in the market.

Q: 2. What is the price of your products?

A: The price of our products is USD18 per unit, with a minimum order of 20 units.

Q: 3. Do your products meet ISO, EN, and CE standards?

A: Yes, our products meet CE, ISO13485, and EN standards. We also have a Nelson certification.

Q: 4. How long does it take to receive the products?

A: It takes 3-7 days to receive the products.

The happy time with customers

