

PE 100% Virgin PE, High Density, Low Density, Linear, Branched, etc.

Product Name:

1. MOQ: 100,000 PCS
2. OEM: Yes
3. Lead Time: 15-20 Days
4. Packaging: 250 PCS/Bag, 51 * 30 * 23 CM/Bag
5. Certificate: FDA, CE, ISO13485
6. Payment: T/T, L/C, etc.

Technical Specification:

Material	PE (High Density, Low Density, Linear, Branched, etc.)
Material	LDPE
Size	50 "x 80" (adult); 40x60"(kids)
Thickness	0.015-0.04 mm.
Color	White
Weight	opp/1 250pcs/bag, 51 * 30 * 23 cm./bag
Quantity	11 bags/lot
MOQ	100,000 PCS
Certificate	FDA, CE, ISO13485

