

Product details

Information

Material: SMS 55g

Weight: 1pc /, 50PCS / CTN

1-3

10,000

1. CE ISO FDA
2. 10,000 / 1 * 40 HQ
3. CE ISO FDA
4. 24

Table

| | |
|----------|----------------------|
| Material | SMS PP, PP / PE, 55g |
| Weight | 1pc /, 50PCS / CTN |
| Material | Snaps |
| Size | S-6XL |
| Material | 5pcs / 50 |
| Material | |
| Material | |

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1.Are you a professional manufacturer of high quality products?

A: We are a professional manufacturer of high quality products since 2003. We have a complete production line and a professional R&D team. We have passed ISO9001:2008 certification in 2013.

Q: 2.What is your minimum order quantity (MOQ) for each product?

A: The MOQ is USD18 per piece. We can accept 20 pieces for each product.

Q: 3. What are the main standards you follow for your products? **ISO, EN** or **CE (s)?**

A: We follow CE, ISO13485, and other international standards. We also have a professional testing center, NELSON.

Q: 4. How do you ensure the quality of your products? Do you have any quality control system?

A: 3-7

The happy time with customers

