

□□□□□ □□□:

□□ □□□□□ □□□ SPP□□□□

□□□□□

1. Dust□□□□□□□□ □□□□□ □□ □□□□□□□□
2. perfect □□□□□ □□□ □□□□□□ □□□□□□□□.
3. □□□□□□□ □□□□ □□□□ □□□□□ □□□□□ □□.
4. Nose□□□□□ □□□□□□□ □ □□□□□□□□ □□□□□□□ □□□ □□ □ □□□□□□ □□□□□□ □□□□□ □□□ □□ □□□□ □□□□.

□□□□□□□□□ □□□□□:

□□□□□□□□□□ □□□□□□□□□□ □ □□□□□□□ □□□□□□ □ □□□□□□□□ □ □□□□□ □□ □□□□□ □□□ □□□□□□□ □□ □□□ □□□□□□ □□□□□

□□□□□□□:

50PCS / 1000 □□□□□□□□□□ / □□□□□□□□

□□□□□□ □□□□ □□□□□□□:

□□□□	SBPP (□□□□□□□□□□□□□□ □□□□□□□□ □□□□)
□□□□□□□ □□□	□□□□
□□□□□□	□□□□ □□□□
□□□□	□□□□□□
□□□□□□	120 □□□□ □□gsm 220 □□□□ gsm
□□□□□□	□□□□□
□□□□□□□□	1000 □ □□□□□ / □□□□□ □□□□□□□□ □□□□ 50 :□□□□□□pcs / ctn □□□□□□ □□□□ □□
□□□□□□	CE□ ISO13485□ FDA
□□□□	100000PCS
□□□□	□□□□ 2 □□ □□□□□□ □□□□□ □□ □□□□□□□□ □□□□
□□□□□□	□□□□□□□□□ □□□ □□□□□□□ □□□□ PP



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INSTRUMENTS

