

The feature of surgical lab coat

-

- ü Isolate external liquids
- ü Filter fine particles and dust
- ü Comfortable and breathable
- ü Suitable for men and women

ü Environmental protection fabric

Fabric Material	SMS
Style	knitted cuffs and collar ,with buttons
Dimensions size	S-3XL
Fabric weight	35-65g, could be as your request
Color	Blue color
Certificate	ISO13485, FDA
MOQ:	10000 pcs
Packing design	can be as your design
Packing	1pc/plastic paper bag,50pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	10-30 days
Applications	laboratory, clinic, food industry, dust-free workshop,hospital

The picture of Surgical Lab coat with knitted cuffs medical supplier




orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL







orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

