

Xiantao supplier for disposable PP heavy weight patient slide sheet CE approved

The specification for heavy weight patient slide sheet

Material	polypropylene
Color	blue
Style	Ultrasonic heat sealing ,with hole in the length side
Size	195cm*85cm (can be made according to client's demand)
Application	Hospital, emergency center, clinic
MOQ	5,000 Pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	10,000pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	According to the quantity, normally 30days after deposit
Packing details	Sterile: 1 Pcs/paper plastic bag,50 Pcs/Ctn Non-sterilize: 1 Pcs/bag, 50 Pcs/ctn Or as your requirements

Features of PP patient slide sheet

- 1) Made of high quality nonwoven fabrics, non woven SMS,PP, or PP+PE 60g-120g;
- 2) With handle, easily used for transferring the patient body;
- 3) Supports up to 300 pounds. Size: 85cm W x 195cm L

Picture of PP heavy weight patient slide sheet



orient Honest
MEDICAL INTERNATIONAL

