

### تفاصيل المنتج:

1. Material: Non woven polypropylene 100%
2. Size:80\*190+20cm,90\*200+20cm,160\*190+20cm,160\*200+20cm,180\*200+20cm or customized size
3. Elastic: 5mm along side width+30cm on length or customized size
4. Colour: Blue; أبيض ؛ other size available
5. Packing: 10pcs/bag, 100pcs/carton

### دواعي الإستعمال:

- 1, Widely used for SPA, Hotel, Massage, Salon, Medical one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

### تحديد:

مادة	Non woven polypropylene 100%
اللون	أبيض أزرق
أسلوب	elastics on width
مقاس	80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
التعبئة	10 pcs/box, 100pcs/ctn
وزن	25-65gsm or customized
طرد	PE bag coated, then paper carton package or customized
ميزة	Durable, soft, disposable
Usage	SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home

### خدماتنا:

1. Provide the highest quality products  
Professional QC team strictly control each step of the production until loading.
2. The fastest delivery time  
If urgent, 1\*40HQ can be finished within 5-7days!  
Normally, 1\*40HQ can be finished within 10-15 days!  
Busy season, 1\*40HQ can be finished within 25-30 days!
3. The most reasonable price  
Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!
- 4.The sweet and quick service  
Our workers will be always so patiently making a response of your any questions within 5-10 minutes.
- 5.We can also help to source connected products.  
Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

