

All isolation gowns are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

Fabric Material	PP+PE material
Style	With ties on collar and waist, knitted cuffs
Dimensions size	X-3XL
Fabric weight	30-40g, could be as your request
Color	Green color
Certificate	ISO13485, FDA
MOQ:	10000 pcs
Packing design	can be as your design
Packing	1pc/plastic paper bag,50pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	7-30 days
<u>Applications</u>	laboratory, clinic, food industry, dust-free workshop,hospital

The picture of disposable Level 1 isolation gown with knitted cuffs medical manufacturer







