

**Produktdetails:**

1. Material: Non woven polypropylene 100%
2. Size:80\*190+20cm,90\*200+20cm,160\*190+20cm,160\*200+20cm,180\*200+20cm or customized size
3. Elastic: 5mm along side width+30cm on length or customized size
4. Colour: Blue; Weiß ; other size available
5. Packing: 10pcs/bag, 100pcs/carton

**Indications**

- 1, Widely used for SPA, Hotel, Massage, Salon, Medical one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

**Spezifikationen:**

Material	Non woven polypropylene 100%
Farbe	Weiß Blau
Stil	elastics on width
Größe	80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
Verpackung	10 pcs/box, 100pcs/ctn
Gewicht	25-65gsm or customized
Paket	PE bag coated, then paper carton package or customized
Merkmal	Durable, soft, disposable
Verwendungszweck	SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home

**Unsere Leistungen:**

1. Provide the highest quality products  
Professional QC team strictly control each step of the production until loading.
2. The fastest delivery time  
If urgent, 1\*40HQ can be finished within 5-7days!  
Normally, 1\*40HQ can be finished within 10-15 days!  
Busy season, 1\*40HQ can be finished within 25-30 days!
3. The most reasonable price  
Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!
4. The sweet and quick service  
Our workers will be always so patiently making a response of your any questions within 5-10 minutes.
5. We can also help to source connected products.  
Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!





