

Nonwoven Clip Cap, Disposable Cap, Mob Cap, Doctor Cap, Bouffant Cap, Nurse Cap, Polypropylene Cap Disposable PP Clip Cap ,Disposable Bouffant Cap

1. Directly manufacturer, disposable clip caps, any colors.
- 2, ISO, CE, good protect, latex free, easy breath, eco friendly,
3. Capacity: 200, 000PCS/day, 1*40HQ can be finished within 10-15 days;
4. OEM packing design, service in 24 hours.

Nonwoven Disposable Caps, Disposable Nurse Cap, Surgical Bouffant Cap With Single Elastic

It is made of PP non woven, sewing and automatic style.

Advantages

1. Good protective, light weight, latex free.
2. Sewing style, Comfortable wearing.
3. White, blue, green, yellow, pink, orange, red colors can be available.

Indications

For hospital clinical medical and patients, dust free work shop, food factory.

Clip caps packing

100PCS/bag, 2000PCS/carton. Carton size: 58X25X27cm 1*20 container can load 715ctns, total 1430, 000PCS. 1*40 container can load 1737ctns, total 3474, 000PCS

| | |
|----------------|---|
| Name | Disposable surgical non woven elastic caps |
| Material | PP 10-30g |
| Style | Bouffant cap, clip cap, mob cap. mop cap. doctor cap |
| Size | 18",19",20",21",22",24" |
| Color | white,green,blue,yellow,pink,red,orange etc. |
| packing | 100 pcs/bag,2000pcs/ctn |
| packing design | all inner bag and carton printing can be as your design |
| MOQ | 100,000 pcs |
| Sample | can be supplied free for your quality checking within 3 days |
| Payment way | T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW, |
| Delivery time | Capacity: 200,000pcs/day, 1*40HQ can be finished within 10-15 days; |
| FOB port | Wuhan/Shanghai/Ningbo/Guangzhou |

The picture of disposable bouffant cap