

## **The feature of surgical lab coat**

- ü Isolate external liquids
- ü Filter fine particles and dust
- ü Comfortable and breathable
- ü Suitable for men and women

### **ü Environmental protection fabric**

Fabric Material	SMS
Style	knitted cuffs and collar ,with buttons
Dimensions size	S-3XL
Fabric weight	35-65g, could be as your request
Color	Blue color
Certificate	ISO13485, FDA
MOQ:	10000 pcs
Packing design	can be as your design
Packing	1pc/plastic paper bag,50pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	10-30 days
Applications	laboratory, clinic, food industry, dust-free workshop,hospital

## **The picture of Surgical Lab coat with knitted cuffs medical supplier**




orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

